



ANG1 Antibody

Product Code	CSB-PA14679A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10152
Immunogen	Recombinant Bovine Angiogenin-1 protein (24-148AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiogenin-1 (EC 3.1.27.-), ANG1, ANG
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	ANG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.