





ASGR1 Antibody

Product Code	CSB-PA17469A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07306
Immunogen	Recombinant Human Asialoglycoprotein receptor 1 protein (68-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Asialoglycoprotein receptor 1 (ASGP-R 1) (ASGPR 1) (C-type lectin domain family 4 member H1) (Hepatic lectin H1) (HL-1), ASGR1, CLEC4H1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ASGR1

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA17469A0Rb at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA17469A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.