



LDHA Antibody, Biotin conjugated

Product Code	CSB-PA00045D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00338
Immunogen	Recombinant Human L-lactate dehydrogenase A chain protein (5-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-lactate dehydrogenase A chain (LDH-A) (EC 1.1.1.27) (Cell proliferation-inducing gene 19 protein) (LDH muscle subunit) (LDH-M) (Renal carcinoma antigen NY-REN-59), LDHA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LDHA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.