





TEP1 Antibody

Product Code	CSB-PA18049A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99973
Immunogen	Recombinant Human Telomerase protein component 1 protein (2-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Telomerase protein component 1 (Telomerase-associated protein 1) (Telomerase protein 1) (p240) (p80 telomerase homolog), TEP1, TLP1, TP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TEP1


Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA18049A0Rb at dilution of 1: 100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA18049A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human esophageal tissue using CSB-PA18049A0Rb at dilution of 1: 100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.