



ANGPT2 Antibody, Biotin conjugated

Product Code	CSB-PA05749D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15123
Immunogen	Recombinant Human Angiopoietin-2 protein (19-485AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-2 (ANG-2), ANGPT2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ANGPT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.