



ILK Antibody, FITC conjugated

Product Code	CSB-PA00999C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13418
Immunogen	Recombinant Human Integrin-linked protein kinase protein (1-228AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin-linked protein kinase (EC 2.7.11.1) (59 kDa serine/threonine-protein kinase) (ILK-1) (ILK-2) (p59ILK), ILK, ILK1 ILK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ILK