





PTMA Antibody

Product Code	CSB-PA01524A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06454
Immunogen	Recombinant Human Prothymosin alpha protein (5-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prothymosin alpha [Cleaved into: Prothymosin alpha, N-terminally processed; Thymosin alpha-1], PTMA, TMSA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTMA

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA01524A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA01524A0Rb at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.