





CSF3 Antibody

Product Code	CSB-PA06459A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09919
Immunogen	Recombinant Human Granulocyte colony-stimulating factor protein (27-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granulocyte colony-stimulating factor (G-CSF) (Pluripoietin) (Filgrastim) (Lenograstim), CSF3, C17orf33 GCSF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CSF3

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA06459A0Rb at dilution of 1:20



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA06459A0Rb at dilution of 1:20



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.