




HLA-DQA2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA14829A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P01906 |
| Immunogen | Recombinant Human HLA class II histocompatibility antigen, DQ alpha 2 chain protein (24-214AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:10-1:100 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | HLA class II histocompatibility antigen, DQ alpha 2 chain (DX alpha chain) (HLA class II histocompatibility antigen, DQ(6) alpha chain) (HLA-DQA1) (MHC class II DQA2), HLA-DQA2, HLA-DXA |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | HLA-DQA2 |

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA14829A0Rb at dilution of 1:20

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.