




# HLA-DQA1 Antibody

<b>Product Code</b>	CSB-PA14839A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P01909
<b>Immunogen</b>	Recombinant Human HLA class II histocompatibility antigen, DQ alpha 1 chain protein (24-213AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	HLA class II histocompatibility antigen, DQ alpha 1 chain (DC-1 alpha chain) (DC-alpha) (HLA-DCA) (MHC class II DQA1), HLA-DQA1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Immunology
<b>Target Names</b>	HLA-DQA1

## Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA14839A0Rb at dilution of 1:50

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.