




eutC Antibody, HRP conjugated

Product Code	CSB-PA18239B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19636
Immunogen	Recombinant Escherichia coli Ethanolamine ammonia-lyase light chain protein (1-295AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ethanolamine ammonia-lyase light chain (EC 4.3.1.7) (Ethanolamine ammonia-lyase small subunit), eutC
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	eutC

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg, 10µg, 5µg

All lanes: eutC antibody; HRP conjugated at 0.42µg/ml

Predicted band size: 32 kDa

Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.