




glpK Antibody, Biotin conjugated

Product Code	CSB-PA18407D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A6F3
Immunogen	Recombinant Escherichia coli Glycerol kinase protein (2-502AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerol kinase (EC 2.7.1.30) (ATP:glycerol 3-phosphotransferase) (Glycerokinase) (GK), glpK
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	glpK

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg, 10µg, 5µg

All lanes: glpK antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 57 kDa

Observed band size: 57 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.