



Mmp13 Antibody, Biotin conjugated

Product Code	CSB-PA17739D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33435
Immunogen	Recombinant Mouse Collagenase 3 protein (105-472AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagenase 3 (EC 3.4.24.-) (Matrix metalloproteinase-13) (MMP-13), Mmp13
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Mmp13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.