



Angpt2 Antibody, Biotin conjugated

Product Code	CSB-PA16959D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35608
Immunogen	Recombinant Mouse Angiopoietin-2 protein (19-483AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-2 (ANG-2), Angpt2, Agpt2
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Angpt2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.