



IGF1 Antibody, HRP conjugated

Product Code	CSB-PA08399B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05019
Immunogen	Recombinant Human Insulin-like growth factor I protein (49-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-like growth factor I (IGF-I) (Mechano growth factor) (MGF) (Somatomedin-C), IGF1, IBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IGF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.