



The metabolite of furazolidone Monoclonal Antibody

Product Code	CSB-MA002511I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	The metabolite of furazolidone-BSA conjugate
Raised In	Mouse
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%?protein G purified
Isotype	IgG1
Clonality	Monoclonal
Alias	CPAOZ
Product Type	Small Molecules Antibodies
Target Names	AOZ
Accession NO.	4G2D7