



HAVCR2 Antibody, Biotin conjugated

Product Code	CSB-PA010145LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDQ0
Immunogen	Recombinant Human Hepatitis A virus cellular receptor 2 protein (22-202AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatitis A virus cellular receptor 2 (HAVcr-2) (T-cell immunoglobulin and mucin domain-containing protein 3) (TIMD-3) (T-cell immunoglobulin mucin receptor 3) (TIM-3) (T-cell membrane protein 3), HAVCR2, TIM3 TIMD3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HAVCR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.