



HAVCR2 Antibody, Biotin conjugated

Product Code	CSB-PA010145LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDQ0
Immunogen	Recombinant Human Hepatitis A virus cellular receptor 2 protein (22-202AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatitis A virus cellular receptor 2 (HAVcr-2) (T-cell immunoglobulin and mucin domain-containing protein 3) (TIMD-3) (T-cell immunoglobulin mucin receptor 3) (TIM-3) (T-cell membrane protein 3), HAVCR2, TIM3 TIMD3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HAVCR2