





DNMT3A Antibody

Product Code	CSB-PA007083LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K1
Immunogen	Recombinant Human DNA (cytosine-5)-methyltransferase 3A protein (680-902AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA (cytosine-5)-methyltransferase 3A (Dnmt3a) (EC 2.1.1.37) (DNA methyltransferase HsaIIIA) (DNA MTase HsaIIIA) (M.HsaIIIA), DNMT3A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DNMT3A

Image



Immunohistochemistry of paraffin-embedded mouse liver tissue using CSB-PA007083LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA007083LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.