




VAMP7 Antibody

Product Code	CSB-PA025785LA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5ZL74
Immunogen	Recombinant Chicken Vesicle-associated membrane protein 7 protein (2-181AA)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle-associated membrane protein 7 (Synaptobrevin-like protein 1), VAMP7, SYBL1
Immunogen Species	Chicken
Research Area	Others
Target Names	VAMP7

Image



Immunohistochemistry of paraffin-embedded mouse liver tissue using CSB-PA025785LA01CH at dilution of 1:50