



GUCA2B Antibody

Product Code	CSB-PA613693GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16661
Immunogen	Human uroguanylin
Raised In	Rabbit
Species Reactivity	human,Mouse,Rat
Tested Applications	ELISA,WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	GUCA2B; GCAP-II; UGN
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GUCA2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.