



BDNF Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA05775D0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23560 |
| Immunogen | Recombinant Human Brain-derived neurotrophic factor protein (137-237AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Brain-derived neurotrophic factor (BDNF) (Abrineurin), BDNF |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | BDNF |