



RAB1A Antibody, Biotin conjugated

Product Code	CSB-PA00335D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62820
Immunogen	Recombinant Human Ras-related protein Rab-1A protein (2-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-1A (YPT1-related protein), RAB1A, RAB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAB1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.