



H3F3A Antibody, Biotin conjugated

Product Code	CSB-PA01435D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84243
Immunogen	Recombinant Human Histone H3.3 protein (2-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H3.3, H3F3A; H3F3B, H3.3A H3F3; H3.3B
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	H3F3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.