



ODC1 Antibody, Biotin conjugated

Product Code	CSB-PA00535D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11926
Immunogen	Recombinant Human Ornithine decarboxylase protein (1-461AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ornithine decarboxylase (ODC) (EC 4.1.1.17), ODC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ODC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.