



ENO1 Antibody, Biotin conjugated

Product Code	CSB-PA02395D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06733
Immunogen	Recombinant Human Alpha-enolase protein (2-428AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (C-myc promoter-binding protein) (Enolase 1) (MBP-1) (MPB-1) (Non-neural enolase) (NNE) (Phosphopyruvate hydratase) (Plasminogen-binding protein), ENO1, ENO1L1 MBPB1 MPB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ENO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.