



RPL27 Antibody, HRP conjugated

Product Code	CSB-PA02345B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61353
Immunogen	Recombinant Human 60S ribosomal protein L27 protein (2-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L27 (Large ribosomal subunit protein eL27), RPL27
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL27
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.