




RPS12 Antibody

Product Code	CSB-PA02525A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25398
Immunogen	Recombinant Human 40S ribosomal protein S12 protein (1-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S12 (Small ribosomal subunit protein eS12), RPS12
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS12

Image



Western blot

All lanes: RPS12 antibody at 2µg/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 15 kDa

Observed band size: 90 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.