



ATP6V1E1 Antibody, FITC conjugated

Product Code	CSB-PA02295C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36543
Immunogen	Recombinant Human V-type proton ATPase subunit E 1 protein (2-224AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	V-type proton ATPase subunit E 1 (V-ATPase subunit E 1) (V-ATPase 31 kDa subunit) (p31) (Vacuolar proton pump subunit E 1), ATP6V1E1, ATP6E ATP6E2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ATP6V1E1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.