



VEGFA Antibody, Biotin conjugated

Product Code	CSB-PA025833LD01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49151
Immunogen	Recombinant Pig Vascular endothelial growth factor A protein (27-190AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor A (VEGF-A) (Vascular permeability factor) (VPF), VEGFA, VEGF
Immunogen Species	Pig
Research Area	Cardiovascular
Target Names	VEGFA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.