




Gstp1 Antibody

Product Code	CSB-PA009989HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04906
Immunogen	Recombinant Rat Glutathione S-transferase P protein (2-210AA)
Raised In	Rabbit
Species Reactivity	Rat, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase P (EC 2.5.1.18) (Chain 7) (GST 7-7) (GST class-pi), Gstp1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Gstp1

Image



Western blot
 All lanes: Gstp1 antibody at 2µg/ml + HeLa whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 24 kDa
 Observed band size: 24 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.