





DEFB1 Antibody

Product Code	CSB-PA09854A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60022
Immunogen	Recombinant Human Beta-defensin 1 protein (33-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 1 (BD-1) (hBD-1) (Defensin, beta 1), DEFB1, BD1 HBD1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	DEFB1


Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA09854A0Rb at dilution of 1:50



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA09854A0Rb at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using CSB-PA09854A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.