





DCN Antibody

Product Code	CSB-PA10289A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07585
Immunogen	Recombinant Human Decorin protein (20-359AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Decorin (Bone proteoglycan II) (PG-S2) (PG40), DCN, SLRR1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DCN

Image



Western blot
 All lanes: DCN antibody at 2µg/ml
 Lane 1: Mouse heart tissue
 Lane 2: HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 40, 28, 24, 20, 9 kDa
 Observed band size: 40 kDa



Immunohistochemistry of paraffin-embedded mouse heart tissue using CSB-PA10289A0Rb at dilution of 1:50



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.