




FOS Antibody

Product Code	CSB-PA008790KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01100
Immunogen	Center-peptide of Human FOS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Storage Buffer	Store at -20oC or -80oC. Avoid freeze / thaw cycles. Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azi
Purification Method	Affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOS

Image



Western blot analysis of extracts of mouse intestine, using c-Fos antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.