




CEACAM5 Antibody

Product Code	CSB-PA005165KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06731
Immunogen	Recombinant protein of Human CEACAM5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:2000, IHC: 1:50-1:200
Relevance	Carcinoembryonic antigen (CEA), also known as CD66e or CEACAM5, is a 180-200 kDa cell surface glycoprotein whose expression is elevated in intestinal carcinomas and other tumors. CEA mediates cell adhesion, though little more is known about its biological activity. Expression of CEA is correlated with tumorigenicity, and it has been shown to play a role in cell migration, adhesion and invasion in culture cells, as well as in metastasis in vivo.
Storage Buffer	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Alias	CEACAM5; CD66e; CEA; DKFZp781M2392
Product Type	Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Intended Use	For research use only. Not for human, diagnostic or therapeutic use.
Target Names	CEACAM5

Image



Western blot analysis of extracts of various cell lines, using CEACAM5 antibody.