




# CYR61 Antibody

<b>Product Code</b>	CSB-PA006463KA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O00622
<b>Immunogen</b>	Recombinant protein of Human CYR61
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; WB: 1:500-1:2000, IHC: 1:50-1:200
<b>Storage Buffer</b>	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	CYR61; IGF1BP1; IGF binding protein 10; IBP 10; CCN family member 1; CCN1
<b>Product Type</b>	Rabbit Anti Human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Intended Use</b>	For research use only. Not for human, diagnostic or therapeutic use.
<b>Target Names</b>	CYR61

## Image



Western blot analysis of Jurkat cell lysate using CYR61 antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.