




CYR61 Antibody

Product Code	CSB-PA006463KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00622
Immunogen	Recombinant protein of Human CYR61
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:2000, IHC: 1:50-1:200
Storage Buffer	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Alias	CYR61; IGFBP10; IGF binding protein 10; IBP 10; CCN family member 1; CCN1
Product Type	Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Intended Use	For research use only. Not for human, diagnostic or therapeutic use.
Target Names	CYR61

Image



Western blot analysis of Jurkat cell lysate using CYR61 antibody.