




# FIGF Antibody

<b>Product Code</b>	CSB-PA008674KA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43915
<b>Immunogen</b>	Recombinant protein of Human FIGF
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
<b>Form</b>	Liquid
<b>Storage Buffer</b>	Store at -20oC or -80oC. Avoid freeze / thaw cycles. Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	Affinity purification
<b>Isotype</b>	IgG
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cancer
<b>Target Names</b>	VEGFD

## Image



Western blot analysis of 293T cell lysate ,using FIGF antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.