




CGA Antibody

Product Code	CSB-PA005293KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01215
Immunogen	Recombinant protein of Human CGA
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Storage Buffer	Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Alias	HCG; LHA; FSHA; GPHa; TSHA; GPHA1; CG-ALPHA
Immunogen Species	Homo sapiens (Human)
Intended Use	For research use only. Not for human, diagnostic or therapeutic use.
Target Names	CGA

Image



Western blot analysis of JAR cell lysate ,using CGA antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.