




IGFBP5 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA011099KA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P24593 |
| Immunogen | Recombinant protein of Human IGFBP5 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200 |
| Storage Buffer | Store at -20°C or -80°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| Purification Method | Affinity purification |
| Isotype | IgG |
| Alias | IGFBP5; IBP5 |
| Product Type | Rabbit Anti Human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Intended Use | For research use only. Not for human, diagnostic or therapeutic use. |
| Target Names | IGFBP5 |

Image



Western blot analysis of extracts of HEK-293T cell lines, using IGFBP5 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.