




IL13 Antibody

Product Code	CSB-PA011590KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35225
Immunogen	Recombinant protein of Human IL13
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Storage Buffer	Store at -20oC or -80oC. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL13

Image



Western blot analysis of extracts of recombinant rotein using IL13 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.