




WNT5A Antibody

Product Code	CSB-PA026138KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41221
Immunogen	Recombinant protein of Human WNT5A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	WNT5A

Image



Western blot analysis of extracts of mouse heart cell lines, using Wnt5a antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.