



PPIA Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA07814C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P62937 |
| Immunogen | Recombinant Human Peptidyl-prolyl cis-trans isomerase A protein (1-165AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peptidyl-prolyl cis-trans isomerase A (PPIase A) (EC 5.2.1.8) (Cyclophilin A) (Cyclosporin A-binding protein) (Rotamase A) [Cleaved into: Peptidyl-prolyl cis-trans isomerase A, N-terminally processed], PPIA, CYPA |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | PPIA |