




PPIA Antibody

Product Code	CSB-PA07814A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62937
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase A protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase A (PPIase A) (EC 5.2.1.8) (Cyclophilin A) (Cyclosporin A-binding protein) (Rotamase A) [Cleaved into: Peptidyl-prolyl cis-trans isomerase A, N-terminally processed], PPIA, CYPA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PPIA

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, Mouse kidney tissue, Mouse brain tissue
 All lanes: PPIA antibody at 4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 19, 12 kDa
 Observed band size: 19 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.