



FN1 Antibody, Biotin conjugated

Product Code	CSB-PA13187D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02751
Immunogen	Recombinant Human Fibronectin protein (627-688AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibronectin (FN) (Cold-insoluble globulin) (CIG) [Cleaved into: Anastellin; Ugl-Y1; Ugl-Y2; Ugl-Y3], FN1, FN
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.