



FOXP3 Antibody, Biotin conjugated

Product Code	CSB-PA15229D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZS1
Immunogen	Recombinant Human Forkhead box protein P3 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein P3 (Scurfin) [Cleaved into: Forkhead box protein P3, C-terminally processed; Forkhead box protein P3 41 kDa form], FOXP3, IPEX
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FOXP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.