



Recombinant Human Hepatocyte nuclear factor 4-alpha (HNF4A)

Product Code	CSB-MP010594HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41235
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRLSKTLVDM DMADYSAALD PAYTTLEFEN VQVLTMGNDT SPSEGTLNLA PNSLGVSALC AICGDRATGK HYGASSCDGC KGFFRRSVRK NHMYSCRFSR QCVVDKDKRN QCRYCRLKKC FRAGMKKEAV QNERDRISTR RSSYEDSSLP SINALLQAEV LSRQITSPVS GINGDIRAKK IASIADVCS MKEQLLVLE WAKYIPAFCE LPLDDQVALL RAHAGEHLLL GATKRSMVFK DVLLLGNDYI VPRHCPELAE MSRVSIRILD ELVLPFQELQ IDDNEYAYLK AIIFDPAK GLSDPGKIKR LRSQVQSLE DYINDRQYDS RGRFGELLLL LPTLQSITWQ MIEQIQFIKL FGMAKIDNLL QEMLLGGSPS DAPHAHHPLH PHLMQEHMGT NVIVANTMPT HLSNGQMCEW PRPRGQAATP ETPQPSPGG SGSEPYKLLP GAVATIVKPL SAIPQPTITK QEVI
Source	Mammalian cell
Target Names	HNF4A
Protein Names	Recommended name: Hepatocyte nuclear factor 4-alpha Short name= HNF-4-alpha Alternative name(s): Nuclear receptor subfamily 2 group A member 1 Transcription factor 14 Short name= TCF-14 Transcription factor HNF-4
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.