



Itgal Antibody

Product Code	CSB-PA011875LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24063
Immunogen	Recombinant Mouse Integrin alpha-L protein (153-325AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-L (CD11 antigen-like family member A) (Leukocyte adhesion glycoprotein LFA-1 alpha chain) (LFA-1A) (Leukocyte function-associated molecule 1 alpha chain) (Lymphocyte antigen 15) (Ly-15) (CD antigen CD11a), Itgal, Lfa-1 Ly-15
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Itgal
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.