



Hpx Antibody, HRP conjugated

Product Code	CSB-PA010723LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20059
Immunogen	Recombinant Rat Hemopexin protein (24-460AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemopexin, Hpx
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Intended Use	Hpx
Target Names	Hpx
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.