



Hpx Antibody, HRP conjugated

Product Code	CSB-PA010723LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20059
Immunogen	Recombinant Rat Hemopexin protein (24-460AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemopexin, Hpx
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Intended Use	Hpx
Target Names	Hpx